

Dry Bond

Dental Drying Agent with Adhesive Booster



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Den-Mat Holdings, LLC
1017 W. Central Avenue
Lompoc, CA 93436 USA

©2018 Den-Mat Holdings, LLC. All rights reserved.
750755100 R:3 06/18CM

FOR PROFESSIONAL USE ONLY

A dental drying agent for use with metal, composites and porcelain surfaces. Non-toxic and non-irritating, Dry Bond should be applied after rinsing Etchant or other cleansing and debriding agents.

CAUTION: Wear protective gloves while using this product.

CAUTION: Wear eye protection while using this product.

CAUTION: This product has not been studied in children, pregnant or breast-feeding women.

DIRECTIONS

Dab directly onto surface to quickly remove oily surface film and moisture. A monomolecular layer of silane bonding agent will remain to increase adhesion. Contact with soft tissue may cause patient discomfort.

Symbols Glossary available at:
www.denmat.com/symbols

