# PERFECTEMP<sup>®</sup> II

 **Instructions For Use**

English .....	1
Français .....	10
Español .....	19
Deutsch .....	28
Nederlands .....	37

## CARTRIDGE CONTENTS

1 Cartridge 50ml (76g) each  
15 Mixing tips

## SYRINGE CONTENTS

2 Syringes 5ml (8g) each  
4 Mixing tips

### *Intra-oral Set*

1:30 min.\*  
\*Approximate

### *Final Set*

6:00 min.

### Temporary Crown & Bridge Material

- Fluorescence
- Automix/Autocure
- Bis-acrylic


 Den-Mat Holdings, LLC  
1017 W. Central Avenue  
Lompoc, CA 93436 USA  
[denmat.com](http://denmat.com)

**dm**  
DenMat

## Perfectemp II

Perfectemp II is a biocompatible Bis-acrylic temporary crown and bridge material available in 50ml (76g) dual-barrel cartridges and 5 ml (8g) dual-barrel syringes. The 4:1 automix system provides thorough mixing with minimal waste. The material will autocure with an initial set of approximately 1:30 minutes intra-orally, providing full polymerization without the use of a light source in 6 minutes. Because Perfectemp II exhibits high strength and abrasion resistance, it is ideal for fabricating single or multi-unit provisionals including long span bridges. When extruded into a pre-impression or a plastic template, the material can be placed directly on to the prepared teeth to form the provisional.

## PHASES

## PRODUCT SPECIFICATIONS

### Type of Material

Perfectemp II is a Temporary Crown & Bridge Class of Material. This Autocure material is based on multifunctional methacrylic esters. Ratio base to catalyst is 4:1.

**Relevant Standards:** None

### Typical Composition

Bis-GMA, ethoxyl	32%
Barium glass	44%
Silica	3%
Polyester	9%
Catalyst, stabilizer, pigments	1%

### Working Properties

Working time	1:00 (minutes)*
Initial set, average	1:30 (minutes)*
Total curing time	≤6:00 (minutes)

### Physical properties

Compressive strength (MPa)	≥150
Transverse strength (MPa) ≥60	

\*This is an average and can vary ±10 seconds.

## PRODUCT SPECIFICATIONS (CONTINUED)

### Physical properties (continued)

Diametral tensile strength (MPa)	31
Barcol hardness	≥20
Water sorption (μg/mm <sup>3</sup> )	≥25

**Safety Data** See SDS Information.

**Storage** 15–25°C/59–77°F

**Shelf-life** 24 mos. from date of manufacture

## DIRECTIONS FOR USE

### 1. Pre-impresion / Clear Plastic Matrix

Prior to the extraction or preparation of the tooth, take an impression using polyvinyl siloxane (such as Clear Bite Registration, DenMat). After the impression sets, trim the interdental gingivae to provide bulk for later finishing of the provisional. If the initial shape and contours of the teeth are not ideal, it is advisable to correct them with flowable or composite. In the case of a missing tooth, you may replace it by securing a denture tooth or forming a composite tooth intra-orally or on a stone model. If using a plastic template, pour a stone model into the impression using a Class IV, high compression stone. Use a vacuum forming machine and clear, rigid tray material such as .020 acrylic to fabricate a plastic template.

## 2. Mixing

Be sure to bleed the cartridge or syringe before applying the mixing tip to ensure an even mix. Cartridge: The Perfectemp II cartridges should be used with an automix dispensing gun specifically designed for crown and bridge material. To ensure that the appropriate dispensing gun is being used, check the side of the gun and verify that it is labeled DS 50 4:1/10:1. Syringe: Place mixing tip onto syringe and extrude. Once extruded through the mixing tip, the material is thoroughly mixed and ready to be placed directly into the pre-impression or plastic template.

- a. Select the appropriate Perfectemp II shade.
- b. Dry the prepared teeth, then place a light application of Vaseline®\* or a water soluble separating medium over prepared teeth and any composite material which will come into contact with the Perfectemp II material.
- c. Remove and discard the cap from the Perfectemp II cartridge/syringe.
- d. Ensure that both cartridge/syringe openings are clear, and bleed a small amount of Perfectemp II before placing the mixing tip on the cartridge/syringe.
- e. Align the mixing tip with the corresponding projections on the cartridge/syringe and turn clockwise, securing the tip into place.

\*Vaseline is a registered trademark of Chesebrough-Pond's Inc.

- f. Squeeze the trigger on the dispensing gun or directly extrude material through dual-barrel syringe into the impression or the clear plastic template.

**Note:** To prevent voids, stir a small instrument inside the mixture and always try to keep the tip in the expressed material as you extrude it moving from tooth to tooth.

- g. Place loaded tray or template within the allotted 1:00 minute\* working time.
- h. Once the template and provisional material has been inserted into the mouth, hold firmly over the prepared teeth. Remove any excess material with a plastic instrument along the gingival margin.
- i. Store the Perfectemp II cartridge/syringe containing the remaining unused material with the used mixing tip intact to avoid any cross contamination of the catalyst and base materials. Do not replace the mixing tip until the next use.

### 3. Removal from the mouth

Remove from mouth at approximately 2:30 minutes\* from start of mix.

## 4. Finishing

Once the provisional has been removed from the mouth, it will become gradually less pliable for an additional 3:30 minutes.\* Allow full polymerization of the material before polishing or contouring.

- a. Bulk trimming, thinning and contouring of the margins and embrasures may be accomplished by using a slow-speed acrylic bur and disk.
- b. A high-speed finishing diamond bur may be used to adjust the occlusion and further correct embrasures while the provisional is in the mouth.  
**Note:** Wear protective eyewear, if necessary, to shield eyes from the composite dust caused during polishing.
- c. An inhibition layer may remain on the surface of the Perfectemp II provisional. This can be removed by wiping with isopropyl alcohol, buffing with a soft fabric wheel, or sand blasting.

## 5. Cementation of the Provisional

Prior to cementing, remove any separating medium on prepared teeth and the provisional. Cement the provisional using resin or zinc oxide - resin based cement.

\*This is an average and can vary  $\pm 10$  seconds.

## ADDING TO, OR REPAIRING PROVISIONALS

These materials are designed to provide a perfect shade match and should be added to a clean and roughened surface after removal of the oxygen inhibited layer. Place flowable composite in desired area, re-seat the restoration and cure.

## STORAGE

Do not store Perfectemp II at temperatures exceeding 25° C / 77° F. Perfectemp II will remain stable for up to 24 months after date of manufacture if kept below 25° C / 77° F.

## DO

- Always store unused material with the hardened mixing tip intact.
- Always ensure that both cartridge/syringe openings are clear of any particles before each use.

## DON'T

- Do not try to force inadvertently set material through the cartridge/syringe openings.
- Do not use excess material outside of the mouth as an indication of intra-oral curing as the temperature will induce faster curing intra-orally.
- Do not replace the cartridge/syringe cap once it has been removed.
- Do not use after the expiration date.
- Do not use this material if the patient exhibits an allergic response to its use or on patients known to have allergic reactions to Bis-acrylic or other methacrylate materials.

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MDSS GmbH