Splash![®]

instant hydrophilicity

MAX

Directions for Use.....	2
Mode d'emploi	11

 Den-Mat Holdings, LLC
1017 W. Central Avenue
Lompoc, CA 93436 USA

R_X ONLY

MD

LOT

EC REP

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Symbols Glossary available at: denmat.com/symbols

Product Description:

Splash!Max Extra Lite Body (Orange) and Lite Body (Green)
ISO 4823; Type 3: Low Consistency

Splash!Max Medium Body (Blue)
ISO 4823; Type 2: Medium Consistency

Splash!Max Heavy Body (Purple)
ISO 4823; Type 1: High Consistency

Range of Indications:

Splash!Max is ideal for dual viscosity impression techniques both one-step and two-step as well as for a variety of dental procedures including, but not limited to, crowns, bridges, inlays, partial and complete dentures, denture repairs, implants and for forming functional peripheries.

Attention:

When using any hemostatic agents, prior to taking an impression, always rinse and then gently rub the tooth with a moist cotton swab to ensure complete removal.

Working & Setting Times

	Min. Working Time	Total Set Time
Splash!Max Regular Set Time 1:1 Cartridge Viscosities	2 min. 15 sec.	4 min. 30 sec.
Splash!Max Half-Time 1:1 Cartridge Viscosities	1 min. 0 sec.	2 min. 15 sec.
Splash!Max Regular Set Time Jumbo 5:1 Cartridge Viscosities	2 min. 15 sec.	4 min. 30 sec.

Directions for Use

Instructions for 1:1 Cartridge (50mL) Viscosities

1. Load the cartridge into the cartridge dispensing gun.
2. Remove cartridge closure and discard.
3. Before applying the mixing tip, extrude impression material until it exits both chambers at the same rate.
4. Attach the mixing tip on the cartridge and squeeze the dispenser to extrude the required amount of material directly onto the prepared tooth or into an impression tray or delivery syringe.
5. Leave the filled mixing tip on the cartridge after use.
6. Before using the cartridge again, remove and discard the used mixing tip. Before re-applying another mixing tip, be sure to extrude material so the product is exiting both chambers at the same rate.
7. Attach a fresh mixing tip and proceed as usual.

Instructions for 5:1 Jumbo Cartridge (360mL)

Viscosities:

1. Disengage and remove the tamper-proof seal on the cartridge by raising the locking lever (fig.1). Discard tamper-proof seal after removal.
2. Insert cartridge into the automix dispensing machine. Start machine and extrude about ¼ inch of the material to ensure even flow of base and catalyst. Discard the unmixed material.
3. Attach the yellow dynamic mixing tip to the end of the cartridge. Fix the mixing tip in place by applying and turning the yellow rotating ring clockwise (fig. 2). Note: use only DenMat dynamic mixing tips and rotating rings for optimum performance.
4. Fill the tray with the desired amount of material by holding the tray in a tilted position and gently pressing against the dynamic mixer. Keep the tip immersed in the material to avoid air bubbles (fig. 3).
5. Do not remove the used mixing tip from the cartridge body as it acts as a seal.
6. Before using again, release the rotating ring (turn counterclockwise) to remove the used tip and check that the outlets in the activating heads are not blocked. Fit a new mixing tip in place and reapply the rotating ring (turn clockwise) and continue with steps 4 – 5.


fig. 1


fig. 2


fig. 3

Splash!Max Technical Specifications

- Cartridge Mixing Technique (XL, L): Use a 1:1/1:2 dispensing cartridge gun, with static mixing tip provided (4.2-12D Diameter).
- Cartridge Mixing Technique (M, H): Use a 1:1/1:2 dispensing cartridge gun, with static mixing tip provided (6.5-11D Diameter), or 5:1 automix unit and yellow dynamic mixing tip with rotating ring for Jumbo Cartridges.
- Storage Conditions: In a dry place at 15 – 25 °C or 59 – 77 °F and relative humidity of 50% or less. Do not expose to direct heat or sun. Refrigeration will not extend shelf life. Always store used cartridges with mixing tips intact to protect remaining material in cartridges from contamination. Re-using the original cap could cause cross-contamination of base and catalyst components and result in premature setting or hardening of the material at the tip. Store cartridges upright.
- Electroplating: Impressions may be copper or silver plated.
- Casting: The impressions can be poured as soon as 30 minutes after removal from the mouth. Surface tension will be reduced if the impression is briefly rinsed out with a surfactant.

Splash!Max Extra Lite Body

Color: Orange
 Recovery from Deformation: 99.8%
 Strain-in-Compression: 5.0%
 Shore A Hardness:..... 43
 Linear Dimensional Change:0.2%

Splash!Max Lite Body

Color: Green
 Recovery from Deformation: 99.8%
 Strain-in-Compression:6.2%
 Shore A Hardness: 41
 Linear Dimensional Change:0.2%

Splash!Max Medium Body

Color: Blue

Cartridge
(50mL)

Jumbo Cartridge
(360mL)

Recovery from Deformation: 99.8% 99.8%
 Strain-in-Compression:3.2% 3.2%
 Shore A Hardness:60 60
 Linear Dimensional Change:0.2% 0.2%

Splash!Max Heavy Body

Color: Purple

Cartridge
(50mL)

Jumbo Cartridge
(360mL)

Recovery from Deformation:99.7% 99.7%
 Strain-in-Compression:3.0% 3.0%
 Shore A Hardness:.....62 62
 Linear Dimensional Change:0.2% 0.2%

Special Notes & Precautions

- Cartridge material must be used with a 1:1/1:2 dispensing cartridge gun system. Jumbo cartridge material must be used in a 5:1 compatible automix unit with dynamic mixing tips and rotating rings.
- The temperature of Splash!Max affects the working time. Ideally it should be kept at room temperature (15 – 25 °C or 59 – 77 °F). Even within those boundaries, changes in temperature will affect the working time: a 2°C (3.6 °F) temperature increase will shorten working time 10-15 seconds; conversely, a 2 °C (3.6 °F) temperature decrease will lengthen working time 10-15 seconds.
- Always “bleed” the cartridge before applying a fresh mixing tip. This will ensure even dispensing of fresh base and catalyst pastes.
- For even more precise application, intraoral tips are available.
- Solutions applied prior to impression taking, such as retraction solution, may impede the setting reaction of the impression. Therefore, material must be cleaned off thoroughly with water and dried.
- Do not leave remains of the impression material in the sulcus.
- The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded. Possible reactions may include temporary numbness or tingling in the soft tissue.

- When taking an impression adjacent to an existing fixed bridge, make sure all undercuts are sufficiently blocked out.
- Some latex gloves or hand lotions could seriously inhibit the setting of vinyl polysiloxane impression materials. All intraoral areas touched with latex gloves should be rinsed and dried prior to impression taking.
- Avoid contaminating Splash!Max impression materials with chemicals that could inhibit the set (e.g. acrylic and methacrylate residues, latex rubber and sulphur compounds).
- Certain hemostatic agents may interfere with the set of Splash!Max impression materials. Care must be taken to wash and dry contaminated areas before applying the syringe material. All intraoral areas touched by latex gloves should be rinsed and dried prior to impression taking.

The information provided by DenMat is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

For professional use only. Federal Law restricts this medical device to sale by or on the order of a dental professional.

Contra-Indications

There are no known contra-indications provided when the product is used as directed. There are no known harmful reactions or secondary effects on patients and/or dental personnel using vinyl polysiloxane impression materials.

Compatibility

Splash!Max is an addition-reaction type silicone impression material and should not be combined with standard or condensation type silicone materials or polyether materials.

For Professional Use Only

Notice: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Splash!Max Extra Lite Body (Orange) et
Lite Body (Vert)
ISO 4823; Type 3 : Faible consistance

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