

CLINICAL STUDY

Rotadent® is safe for use on implant surfaces

Authors: Thomson-Neal D, Evans GH, Meffert RM

Publication: The International Journal of Periodontics and Restorative Dentistry

Date: 1989

Reference: Thomson-Neal D, Evans GH, Meffert RM. Effects of various prophylactic treatments on titanium, sapphire and hydroxyapatite-coated implants: an SEM study. Intl J Periodontics Restorative Dent. 1989; 9(4): 300 – 311.

Abstract:

The purpose of this study was to evaluate the effects of available professional and personal prophylactic modalities on different implant surfaces.

Methodology:

- In-vitro short-term study
- Four types of dental implants were obtained: two types of titanium implants, hydroxyapatite-coated titanium implants, single crystal sapphire implants

Results:

- The modalities that showed slight or no change included Rotadent on all implant surfaces tested

Conclusion:

Rotadent is a safe home oral hygiene regimen to use on all implant surfaces studied.

