

## Perio Restore

### Homecare Treatment for Periodontal Disease

#### What is Perio Restore?

A homecare treatment system for your patients to keep periodontal disease at bay in-between each dental visit.

The Perio Restore System consists of a customized tray designed to gently push a 1.7% hydrogen peroxide cleansing gel deep into each periodontal pocket, allowing deeper reach than other non-invasive treatments.

- ✓ Each **customized tray** is designed to treat below the gumline and interproximally
- ✓ The oral cleansing gel significantly reduces bacteria with each 10-15 minute application,\* while providing **fresh breath and brighter teeth**
- ✓ **Fights bacteria** inside the pocket and keeps it from recolonizing



“Perio Restore homecare has become a critical extension of what we do in the chair. It’s an effective, non-invasive solution for prevention and treatment of periodontal disease for all of our patients.”

Suzanne Mericle, DMD


\*Eckhard Alt, Francesca Leipold, Danica Milatovic, Günter Lehmann, Sybille Heinz, Albert Schömig, Hydrogen peroxide for prevention of bacterial growth on polymer biomaterials, The Annals of Thoracic Surgery, Volume 68, Issue 6, 1999, Pages 2123-2128, ISSN 0003-4975, [https://doi.org/10.1016/S0003-4975\(99\)00832-2](https://doi.org/10.1016/S0003-4975(99)00832-2). (<http://www.sciencedirect.com/science/article/pii/S0003497599008322>)

## How does Perio Restore treat biofilm?

Periodontal disease is caused by anaerobic bacteria forming plaque in the periodontal pockets which colonizes and forms biofilm.

When perio disease is not controlled, it can lead to other systemic issues such as heart disease, adverse pregnancy outcomes, and diabetes, amongst other issues. An at-home adjunctive solution is needed to treat the disease.

**Solution:** In conjunction with scaling and root planing, the use of Perio Restore is extremely effective at treating periodontal disease. The 1.7% Hydrogen Peroxide gel breaks up biofilm, kills bacteria, and keeps the bacteria from recolonizing, preventing further infection.



## Getting started is simple!

1. Take your patient's digital or VPS impressions
2. Fill out the Perio Restore Lab Rx Form
3. Send Rx Form and impressions digitally or by mail to the DenMat Lab
4. Receive patient case within 9 – 12 business days
5. Send patient home with custom trays and Perio Restore Oral Cleansing Gel (sold separately)



To view the latest news, education and more, visit:

[www.clinicianschoice.com](http://www.clinicianschoice.com)    

© 2025 Clinician's Choice® Dental Products, Inc. All rights reserved.