



Splash!

VPS Impression Material

Splash!® is the impression material for the clinician who demands reliability and value in a VPS impression material. Splash! features excellent flow properties, exceptional tear strength and a wild berry scent that patients will love.

Available in two set times, five high-contrast viscosities and four delivery methods, Splash! is the superhydrophilic impression material for every technique and every impression.

- Super hydrophilic
- Two set times (2:15 and 4:30)
- Available in five high-contrast colors and four delivery methods
- Wild berry scent



| Features | Benefits |
|----------------------------|---|
| Low contact angle | Increased hydrophilicity means more accurate castings |
| High tear strength | Will not tear upon removal from gingival sulcus |
| High dimensional stability | Multiple pours and exceptional detail reproduction indefinitely |
| Wild berry scent | Reduced gag reflex |

| Product | Item # |
|---|---------|
| Splash! 2-Pack Cartridges Regular Set (4:30) | |
| Extra Lite Body | SPD1208 |
| Lite Body | SPD1209 |
| Medium Body | SPD1210 |
| Heavy Body | SPD1211 |
| <i>Contents: 2 – 48 ml cartridges, tips, instructions</i> | |
| Splash! 2-Pack Cartridges Half-Time Set (2:15) | |
| Extra Lite Body | SPD1230 |
| Lite Body | SPD1231 |
| Medium Body | SPD1232 |
| Heavy Body | SPD1233 |
| <i>Contents: 2 – 48 ml cartridges, tips, instructions</i> | |

| Product | Item # |
|---|---------|
| Splash! Bulk Packs Regular Set (4:30) | |
| Extra Lite Body (8-pack) | SPD1587 |
| Lite Body (8-pack) | SPD1588 |
| Medium Body (20-pack) | SPD1589 |
| Heavy Body (20-pack) | SPD1590 |
| <i>Contents: Extra Lite & Lite: 8 – 48 ml cartridges; Medium & Heavy: 20 – 48 ml cartridges</i> | |
| Splash! Bulk Packs Half-Time Set (2:15) | |
| Extra Lite Body (8-pack) | SPD1591 |
| Lite Body (8-pack) | SPD1592 |
| Medium Body (20-pack) | SPD1593 |
| Heavy Body (20-pack) | SPD1594 |
| <i>Contents: Extra Lite & Lite: 8 – 48 ml cartridges; Medium & Heavy: 20 – 48 ml cartridges</i> | |

| Product | Item # |
|--|---------|
| Splash! Jumbo Cartridge Refill Kit Regular Set (4:30) | |
| Heavy Body | SPD1559 |
| <i>Contents: 2 – 380 ml jumbo cartridges, tips, instructions</i> | |
| Splash! Putty Jars | |
| Regular Set Putty Jars (5:30) | SPD1212 |
| Half-Time Set Putty Jars (2:45) | SPD1235 |
| <i>Contents: 1 – 250 ml jar base, 1 – 250 ml jar catalyst, mixing spoons, instructions</i> | |
| Splash! Putty Paks | |
| Regular Set Putty Paks (5:30) | SPD1219 |
| Half-Time Set Putty Paks (2:45) | SPD1237 |
| <i>Contents: 12 – 18 ml paks, instructions</i> | |

To view the latest news, education and more, visit:

www.clinicianschoice.com

