

# SAFETY DATA SHEET

| Issuing Date   | 22-May-2015            | <b>Revision Date</b> | 22-May-2015          | ECR# 23573                  |                |
|--|------------------------|----------------------|----------------------|-----------------------------|----------------|
| 1. IDENT   | <b>IFICATION OF TH</b> | HE SUBSTANC          | CE/PREPARATION       | I AND THE COMPANY/U         | INDERTAKING    |
| GHS product  | <u>identifier</u>      |                      |                      |                             |                |
| Product Name   | 9                      | Solvenex             |                      |                             |                |
| Other means  | of identification      |                      |                      |                             |                |
| Synonyms   |                        | None                 |                      |                             |                |
| Recommende   | ed use of the chemica  | Il and restrictions  | s on use             |                             |                |
| Recommende   | ed Use                 | Professional de      | ental cement remover |                             |                |
| Uses advised   | against                | No information       | available            |                             |                |
| Supplier's det   | ails                   |                      |                      |                             |                |
| Supplier A<br>DenMat<br>1017 W. C<br>Lompoc, C<br>TEL: 805-3 | entral Ave.<br>A 93436 |                      |                      |                             |                |
| Emergency te   | lephone number         |                      |                      |                             |                |
| Emergency Te<br>Number                                       | elephone               | 805-346-3700         |                      |                             |                |
|  |                        | 2. HAZ               | ARDS IDENTIFICA      | ATION                       |                |
| Classificatio  |                        | us according to the  | e OSHA Hazard Comm   | unication Standard 2012 (29 | CFR 1910.1200) |
| Serious Eye Da   | amage/Eye Irritation   |                      |                      | Category 2                  |                |
| GHS Label ele  | ements, including pre  | ecautionary state    | ments_               |                             |                |

**Emergency Overview** 

Warning

Signal Word Hazard Statements • Causes serious eye irritation



#### **Precautionary Statements**

## Prevention

- Wash face, hands and any exposed skin thoroughly after handling.
- Wear eye/face protection.

#### **General Advice**

None

#### Eyes

• IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

• If eye irritation persists: Get medical advice/attention.

### Storage

None

## Disposal

None

## Hazard Not Otherwise Classified (HNOC)

Not applicable

# Other information

May cause skin irritation

# 3. COMPOSITION/INFORMATION ON INGREDIENTS

| Chemical Name   | CAS-No  | Weight % | Trade secret |  |  |
|---|---------|----------|--------------|--|--|
| Citric acid   | 77-92-9 | 30-35    | *            |  |  |
| Isopropyl alcohol   | 67-63-0 | 5-10     | *            |  |  |
| *The superior of a superior (i.e., ) of a superior back have with ball as a fundament ( |         |          |              |  |  |

\*The exact percentage (concentration) of composition has been withheld as a trade secret.

# 4. FIRST AID MEASURES

| Description of necessary first-aid measures |   |  |  |  |
|---|---|--|--|--|
| Eye Contact                                 | Immediately flush with plenty of water. After initial flushing, remove any contact lenses and continue flushing for at least 15 minutes. Get medical attention if irritation persists.  |  |  |  |
| Skin Contact                                | Wash off immediately with soap and plenty of water. Get medical attention if irritation develops and persists.  |  |  |  |
| Inhalation                                  | Move to fresh air. If symptoms persist, call a physician.   |  |  |  |
| Ingestion                                   | Do NOT induce vomiting. Clean mouth with water and afterwards drink plenty of water.<br>Never give anything by mouth to an unconscious person. Get medical attention if symptoms occur. |  |  |  |

# Most important symptoms/effects, acute and delayed

Most Important Symptoms/Effects Red eyes. Tearing. Irritation.

#### Indication of immediate medical attention and special treatment needed, if necessary

Notes to Physician

Treat symptomatically.

# **5. FIRE-FIGHTING MEASURES**

#### Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media No information available.

#### Specific Hazards Arising from the Chemical

Thermal decomposition can lead to release of irritating gases and vapors.

| Explosion Data                   |       |
|----------------------------------|-------|
| Sensitivity to Mechanical Impact | None. |
| Sensitivity to Static Discharge  | None. |

# Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

# 6. ACCIDENTAL RELEASE MEASURES

#### Personal precautions, protective equipment and emergency procedures

| Personal Precautions   | Use personal protective equipment. Ensure adequate ventilation. Keep people away from and upwind of spill/leak. Avoid contact with skin, eyes and clothing. Wash thoroughly after handling. |  |  |
|--|---|--|--|
| Environmental Precautions                                    |   |  |  |
| Environmental Precautions                                    | Avoid release to the environment. Dispose of contents/container to an approved waste disposal plant. See Section 12 for additional Ecological Information.                                  |  |  |
| Methods and materials for contain                            | ment and cleaning up  |  |  |
|  |   |  |  |
| Methods for Containment                                      | Prevent further leakage or spillage if safe to do so.   |  |  |
| Methods for Cleaning Up                                      | Soak up with inert absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly.  |  |  |
|  | 7. HANDLING AND STORAGE   |  |  |
| Precautions for safe handling                                |   |  |  |
| Handling   | Wear personal protective equipment. Ensure adequate ventilation. Do not get in eyes. Avoid contact with skin, eyes and clothing.  |  |  |
| Conditions for safe storage, including any incompatibilities |   |  |  |
| Storage  | Keep containers tightly closed in a dry, cool and well-ventilated place. Keep in properly labeled containers. Keep out of the reach of children.  |  |  |

**Incompatible Products** None known based on information supplied.

# 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

# **Control parameters**

# **Exposure Guidelines**

| Chemical Name                   | ACGIH TLV                     | OSHA PEL  | NIOSH IDLH  |
|---------------------------------|-------------------------------|---|---|
| Isopropyl alcohol<br>67-63-0    | STEL: 400 ppm<br>TWA: 200 ppm | TWA: 400 ppm<br>TWA: 980 mg/m <sup>3</sup><br>(vacated) TWA: 400 ppm<br>(vacated) TWA: 980 mg/m <sup>3</sup><br>(vacated) STEL: 500 ppm | IDLH: 2000 ppm 10% LEL<br>TWA: 980 mg/m <sup>3</sup><br>TWA: 400 ppm<br>STEL: 500 ppm<br>STEL: 1225 mg/m <sup>3</sup> |
| ppropriate engineering controls |                               | (vacated) STEL: 1225 mg/m <sup>3</sup>  | 0   |

| Engineering Measures | Showers             |
|----------------------|---------------------|
|                      | Eyewash stations    |
|                      | Ventilation systems |

# Individual protection measures, such as personal protective equipment

| Eye/Face Protection                                | Safety glasses with side-shields. If splashes are likely to occur, wear: Tightly fitting safety goggles.   |
|--|--|
| Skin and Body Protection<br>Respiratory Protection | Wear protective gloves/clothing.<br>If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved<br>respiratory protection should be worn. Positive-pressure supplied air respirators may be<br>required for high airborne contaminant concentrations. Respiratory protection must be<br>provided in accordance with current local regulations. |
| Hygiene Measures                                   | Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or smoke when using this product. Keep away from food, drink and animal feeding stuffs. Wash hands before breaks and immediately after handling the product. Provide regular cleaning of equipment, work area and clothing.   |

# 9. PHYSICAL AND CHEMICAL PROPERTIES

# Information on basic physical and chemical properties

| Physical State   | Liquid   | Appearance     | Blue                     |
|--|--|----------------|--------------------------|
| Odor   | Mild, Alcohol  | Odor Threshold | No information available |
| Odor<br>Property<br>pH<br>Melting Point/Range<br>Boiling Point/Boiling Range<br>Flash Point<br>Evaporation rate<br>Flammability (solid, gas)<br>Flammability Limits in Air<br>upper flammability limit<br>lower flammability limit<br>lower flammability limit<br>Vapor Pressure<br>Vapor Density<br>Specific Gravity<br>Water Solubility<br>Solubility in other solvents<br>Partition coefficient: n-octand<br>Autoignition Temperature | Values< 3.0No data availableNo data availableInsoluble |                | No information available |
| Decomposition Temperature  | No data available  | None known     |                          |
| Viscosity  | Similar to water   | None known     |                          |
| Flammable Properties   | Not flammable  |                |                          |

#### Explosive Properties Oxidizing Properties

No data available No data available

# Other information

VOC Content (%)

8.67%

# **10. STABILITY AND REACTIVITY**

## **Reactivity**

No data available.

# **Chemical stability**

Stable under recommended storage conditions.

# Possibility of hazardous reactions

None under normal processing.

## Hazardous Polymerization

Hazardous polymerization does not occur.

#### Conditions to avoid

None known based on information supplied.

# Incompatible materials

None known based on information supplied.

# Hazardous decomposition products

Carbon oxides.

# **11. TOXICOLOGICAL INFORMATION**

# Information on likely routes of exposure

| Product Information |  |
|---------------------|--|
| Inhalation          | There is no data available for this product.   |
| Eye Contact         | There is no data available for this product. Expected to be an irritant based on components  |
| Skin Contact        | There is no data available for this product. Prolonged or repeated contact may cause skin<br>irritation  |
| Ingestion           | There is no data available for this product. Ingestion may cause irritation to mucous membranes. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea. |

| Symptoms related to the physical, chemical and toxicological characteristics             |  |  |  |
|--|--|--|--|
| Symptoms   | Eye contact with liquid may cause irritation including stinging, burning, tearing, or reddening of the eyes. |  |  |
| Delayed and immediate effects and also chronic effects from short and long term exposure |  |  |  |
| Sensitization  | No information available.  |  |  |

#### Mutagenic Effects Carcinogenicity

No information available. The table below indicates whether each agency has listed any ingredient as a carcinogen.

| Chemical Name     | ACGIH | IARC    | NTP | OSHA |
|-------------------|-------|---------|-----|------|
| Isopropyl alcohol |       | Group 3 |     |      |

#### IARC: (International Agency for Research on Cancer)

Group 3 - Not Classifiable as to its Carcinogenicity to Humans

| Reproductive Toxicity    | No information available. |
|--------------------------|---------------------------|
| STOT - single exposure   | No information available. |
| STOT - repeated exposure | No information available. |
| Aspiration Hazard        | No information available. |

#### Numerical measures of toxicity - Product

The following values are calculated based on chapter 3.1 of the GHS document:LD50 Oral6253 mg/kg; Acute toxicity estimateLD50 Dermal21376 mg/kg; Acute toxicity estimate

# **12. ECOLOGICAL INFORMATION**

#### Ecotoxicity

The environmental impact of this product has not been fully investigated.

| Chemical Name                | Toxicity to Algae  | Toxicity to Fish   | Toxicity to<br>Microorganisms | Daphnia Magna (Water<br>Flea)              |
|------------------------------|--|--|-------------------------------|--|
| Citric acid<br>77-92-9       |  | LC50 96 h: = 1516 mg/L<br>static (Lepomis macrochirus)   |                               | EC50 72 h: = 120 mg/L<br>(Daphnia magna)   |
| Isopropyl alcohol<br>67-63-0 | EC50 96 h: > 1000 mg/L<br>(Desmodesmus subspicatus)<br>EC50 72 h: > 1000 mg/L<br>(Desmodesmus subspicatus) | LC50 96 h: = 11130 mg/L<br>static (Pimephales promelas)<br>LC50 96 h: = 9640 mg/L<br>flow-through (Pimephales<br>promelas)<br>LC50 96 h: > 1400000 µg/L<br>(Lepomis macrochirus) |                               | EC50 48 h: = 13299 mg/L<br>(Daphnia magna) |

Persistence and Degradability

No information available.

**Bioaccumulation** 

No information available.

| Chemical Name     | Log Pow |  |  |  |
|-------------------|---------|--|--|--|
| Citric acid       | -1.72   |  |  |  |
| Isopropyl alcohol | 0.05    |  |  |  |
|                   |         |  |  |  |

Other Adverse Effects

No information available.

# **13. DISPOSAL CONSIDERATIONS**

| Waste Disposal Methods   | This material, as supplied, is not a hazardous waste according to Federal regulations (40 CFR 261). This material could become a hazardous waste if it is mixed with or otherwise comes in contact with a hazardous waste, if chemical additions are made to this material, or if the material is processed or otherwise altered. Consult 40 CFR 261 to determine whether the altered material is a hazardous waste. Consult the appropriate state, regional, or local regulations for additional requirements. |  |  |  |
|--------------------------|---|--|--|--|
| Contaminated Packaging   | Do not re-use empty containers.   |  |  |  |
| 14 TRANSPORT INFORMATION |   |  |  |  |

DOT

Not regulated

| TDG        | Not regulated |  |  |  |
|------------|---------------|--|--|--|
| <u>MEX</u> | Not regulated |  |  |  |
| ICAO       | Not regulated |  |  |  |
| IATA       | Not regulated |  |  |  |
| IMDG/IMO   | Not regulated |  |  |  |
| <u>RID</u> | Not regulated |  |  |  |
| ADR        | Not regulated |  |  |  |
| ADN        | Not regulated |  |  |  |
|            |               |  |  |  |

# **15. REGULATORY INFORMATION**

| International Inventories |        |
|---------------------------|--------|
| TSCA                      | Exempt |
| DSL                       | Exempt |

#### Legend

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

# U.S. Federal Regulations

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

#### SARA 311/312 Hazard Categories

| Acute Health Hazard               | Yes |
|-----------------------------------|-----|
| Chronic Health Hazard             | No  |
| Fire Hazard                       | No  |
| Sudden Release of Pressure Hazard | No  |
| Reactive Hazard                   | No  |

#### Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

## **CERCLA**

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302):

| Chemical Name          | Hazardous Substances RQs | Extremely Hazardous Substances<br>RQs | RQ                        |
|------------------------|--------------------------|---------------------------------------|---------------------------|
| Citric acid            | -                        | -                                     | RQ Section number 180.950 |
| U.S. State Regulations |                          |                                       |                           |

# California Proposition 65

This product does not contain any Proposition 65 chemicals.

# U.S. State Right-to-Know Regulations

| Chemical Name     | New Jersey | Massachusetts | Pennsylvania | Illinois | Rhode Island |
|-------------------|------------|---------------|--------------|----------|--------------|
| Isopropyl alcohol | Х          | Х             | Х            |          | Х            |

#### U.S. EPA Label Information

# EPA Pesticide Registration Number Not applicable

| 16. OTHER INFORMATION          |  |              |   |                   |                                    |
|--------------------------------|--|--------------|---|-------------------|------------------------------------|
| NFPA                           | Health Hazard 2  | Flammability | 0 | Instability 0     | Physical and Chemical<br>Hazards - |
| HMIS                           | Health Hazard 2  | Flammability | 0 | Physical Hazard 0 | Personal Protection X              |
|                                |  |              |   |                   |                                    |
| Prepared By                    | Product Stewardship<br>23 British American Blvd.<br>Latham, NY 12110<br>1-800-572-6501 |              |   |                   |                                    |
| Issuing Date                   | 22-May-2015  |              |   |                   |                                    |
| Revision Date<br>Revision Note | 22-May-2<br>Initial Rel  |              |   |                   |                                    |

<u>General Disclaimer</u> The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

**End of Safety Data Sheet**